

FACE MEASUREMENTS

FACE MASK & CHIN STRAP – 4A



**Therapist
Support
Laboratory**

PO Box 1022
Collingwood, Victoria, 3066

Tel: 03 9421 2880

Fax: 03 9421 1335

Website: www.tsลาustralia.com.au



Email: tsl@tsลาustralia.com.au

Patient Name _____ **M** | Date _____

OT/Physio _____ **F** | Diagnosis _____

Age _____

Garment Descriptions	Qty	Item No
Chin strap		4001
Chin strap variation		4002
Face mask		4003
Open face mask		4004
Open head face mask		4005
Open head & open face mask		4006
Garment Specifications	Qty	Item No.
Nose covering		4007
Ear covering L/R		4008
Lip covering		4009
Lining (indicate position)		6001
Detached crew neck 1.3 cm		6020
Detached crew neck 2.5cm		6021
Detached crew neck 5.0		6022
Silipos Lining		6024
Silon Tex Lining		60241



A	Circumference at Chin Angle	cm
B	Length of Nose	cm
C	Distance between Eye & Ear	cm
D	Mouth to Chin	cm
E	Tip of Chin to Throat	cm
F	Neck to end of garment	cm
G	Neck Circumference	cm
H	Above Eyebrow to req. height (If open head)	cm
I	Circ. around head at forehead	cm
J	Length of ear	L R cm cm
K	Distance between eyes	cm
L	Length of one eye	cm
M	Under nose to mouth	cm
N	Width under nose	cm
O	Eyebrow To chin	cm



Important points!

Place all measurements in the boxes provided at right.
Circumferential measurements must be taken with a dressmaking measuring tape.
The remaining vertical and horizontal measurements must be taken with a ruler.
All face masks have a zipper fastening unless otherwise indicated.
All chin straps have Velcro fastening.
If lining is required it is recommended you also indicate the position on the facial diagram.